



LESSONS FROM NEANDERTAL MICROBIOTA: WORKING TOWARD ACCURATE METAGENOMIC ANALYSIS FROM LOW-BIOMASS SAMPLES.

LAURA S. WEYRICH
AUSTRALIAN CENTRE FOR ANCIENT DNA
UNIVERSITY OF ADELAIDE, AUSTRALIA




PHOTOGRAPH BY JOE MCNALLY, NATIONAL GEOGRAPHIC

EVOLUTIONARY HISTORY OF THE HOMINID MICROBIOTA




Li et al. BMC Microbiology 2013




Moeller, et al. PNAS, 2014.


ANCIENT DENTAL CALCULUS IS A FOSSILISED BACTERIAL RECORD



Adler et al. Nat Gen, 2013.



Warinner et al. Nat Gen, 2014.



Adler et al. Nat Gen, 2013.

WHY ARE ANCIENT SAMPLES PROBLEMATIC?



PCR

10^6 DNA copies
per
droplet (~ 0.005 uL)




Ancient Sample

$<10^6$ DNA copies
per
1 g of sample

vs

**Low 'endogenous'
DNA yield**

Fragmented and damaged DNA
crosslinking/hydrolysis/oxidation



Willerslev and Cooper. Proc Biol Sci. 2005 Jan 7; 272(1558): 3–16



Contamination
(modern or environmental)
sample collection/handling &
lab or reagent contamination


Modern samples can also be similarly problematic



Skin swabs
Placenta and Pre-term infants
Respiratory biopsies
Blue Ice
Ancient/Historic Specimens
and more...



CONCEPTS IN CONTAMINATION



1. Dedicated DNA extraction/processing facilities
2. Decontamination, if possible
3. Monitoring laboratory & environmental contamination
4. Efficient DNA extractions and library preparation
5. Accurate analysis and verification tools

1. REDUCING CONTAMINATION WITH CLEAN FACILITIES




Isolated pre-PCR facility
Isolated ventilation
Glove boxes
Nightly UV irradiation
Clean -> Dirty workflow
Bleach treatment of stock
Reagents tested for DNA
Appropriate dress

2. DECONTAMINATION MATTERS!

Decontamination Procedures



- A. Remove surface
- B. Decontaminate sample
 - UV or washing, if able
- C. Include controls every time
 - PCR negatives
 - Extraction blank controls
- D. Duplicate samples
- E. Remove contaminants with bioinformatic methods

DNA Contamination




A. Farrer, in prep.


Check out the kitome: Salter, *et al.*, 2014, BMC Biology, 12, 87.



3. DNA EXTRACTION BIASES AND EFFICIENCY MUST BE CONSIDERED



Mackenzie, et al. Front Microbiol. 2015; 6: 130.




Extraction methods bias can results, and limit downstream comparisons.

QIITA!

Extraction efficiency matters!


3. LIBRARY PREPARATION ALSO MATTERS



16S rRNA gene seq will likely not work for highly degraded samples.

Filtering helps (sometimes!)

4. BIOINFORMATIC METHODS TO IDENTIFY SPECIES FROM SAMPLES



Weyrich, et al, Nature, 2017.

Satisfies input requirements; accurate; rapid
nucleotide vs protein

Current methods applied: DIAMOND, MALT

4. BIOINFORMATIC METHODS TO ACCOUNT FOR CONTAMINATION

STRINGENT

1. Decontam <https://github.com/benjjneb/decontam>
2. MEGAN6CE <https://github.com/husonlab/megan-ce>
3. Direct Filtering QIIME2, etc.
4. Contaminant Assessment
Comparison to known lists in:
Salter, et al., *BMC Biol*, 2014 or
Weyrich, et al., *MER*, 2019

RIDE CHECKLIST FOR PERFORMING/REVIEWING LOW MICROBIAL BIOMASS MICROBIOME STUDIES

Report the experimental design and approaches used to reduce and assess the contributions of contamination.

Include controls to assess contaminant DNA. One of each type of negative control (sampling blanks, DNA extraction blanks, and no-template amplification) must be included per sampling, extraction, or amplification batch.

Determine the level of contamination by comparing biological samples to controls.

Explore contaminant taxa within each study and report their impact on the interpretation of biological samples.

WHAT FACTORS SHOULD BE CONSIDERED BEFORE WORKING WITH LOW BIOMASS SAMPLES?

Before you start...

Consider your working environment.

Assess contamination.

Decide on the best extraction method.

During analysis...

Use clean environments for processing.

Include extraction blank and PCR controls.

Don't let bioinformatics be a black box.

Scrutinise and test your results!

And enjoy!